

BLANK PAGE






Indian Standard

SPECIFICATION FOR TREPHINE, BONE LACRIMAL WITH HANDLE

- 1. Scope Covers requirements for four sizes of trephine with handle (with and without central drill point) used for bones lacrimal.
- 2. Dimensions As shown in Fig. 1.



D = 7, 8, 9 OR 10 $D_1 = 5.8, 6.8, 7.8 \text{ OR } 8.8$

All dimensions in millimetres.

FIG. 1 TREPHINE BONE, LACRIMAL (WITH AND WITHOUT CENTRAL DRILL)

- 2.1 A deviation of ± 2.5 percent shall be allowed on all dimensions.
- 3. Material
 - Crown Stainless steel conforming to designation 40Cr13 of IS: 6603-1972 'Specification for stainless steel bars and flats.'
 - Shaft Stainless steel conforming to designation 20Cr13 of IS: 6603-1972.
 - Handle Brass having minimum 60 percent copper.
- **4.** Workmanship and Finish The trephine shall be free from burrs, pits, cracks and other surface defects. The cutting serrations of the working end shall be well formed and uniform. Joint of the shaft and crown shall be neat and sound and done by brazing. The shaft of the trephine shall fit nicely without slippage to the handle. The handle shall be integrated with the trephine.
- 4.1 The handle shall be properly made and shall be free from surface defects. All edges of the handle shall be even and rounded off. The handle shall be plated chromium over nickel and the plating shall conform to Grade 2 of IS: 4827-1968 'Specification for electroplated coatings of nickel and chromium on copper and copper alloys'.
- 4.2 The shape of the handle may be octagonal or hexagonal or as shown in Fig. 1.
- 5. Hardness 600 to 700 HV for crown 400 to 450 HV for shaft

Adopted 28 May 1985 © October 1985, ISI Gr 1

IS: 11244 - 1985

6. Tests

- **6.1** Corrosion Resistance Test The crown and shaft of trephine shall be tested in accordance with IS: 7531-1975 'Method for boiling and autoclaving test for corrosion resistance of stainless steel surgical instruments and shall not show any sign of corrosion after the test.
- **6.2** Performance Test A piece of fresh bone shall be cut with the trephine. The trephine shall cut the bone neatly without getting damaged and shall continue to function in the normal way after the test.
- 7. Marking Each trephine shall be legibly and indelibly marked with the manufacturer's name, initials or recognized trade-mark.
- 7.1 ISI Certification Marking Details available with the Indian Standards Institution.
- 8. Packing As agreed to between the purchaser and the supplier.